Pro-Pul-Pan

Zinc Oxide-Eugenol Cement

Temporary restorative cement

Composition

Poudre Zinc Oxide 90.5 %, Polymethylmethacrylate 5 %, Zinc Acetate 0.5 %, excipient

Liquide Eugenol

Indications

- ▶ Temporary fillings
- ▶ Temporary sealing of crowns and inlays
- ▶ Indirect pulp capping

Directions for use

- ▶ Place powder and liquid on a mixing pad at a recommended powder/liquid (P/L) weight ratio between 1.8 and 3 (g/g)
- ▶ Three spatulation cycles are required to obtain a smooth, creamy mass of the desired consistency
 - Half the powder should be mixed in the first cycle and the remainder in the two subsequent cycles
 - Speed is not essential as setting will only begin 3 or 4 minutes after initiating the mix
 - Mixing on a cool, dry slab induces a slower set. A thin mix also means a slow set
 - Mixing on a warm, damp slab induces a faster set. A thick mix also sets faster
 - The setting time (tested according to ISO 3107) is 10-15 minutes, depending on the P/L ratio.

Storage

Store at room temperature between 5°C and 30°C, protected from direct light.

Warning

Due to the presence of eugenol, avoid using cement with resin-based fillings. When using this product, gloves, mask and protective glasses should be worn. In case of contact with the skin wash with plenty of soap and water. If skin irritation occurs, get medical advice/attention.

In case of contact with the eye rinse cautiously with soft, clean water for 15 minutes, holding the eyelids open, and call a doctor immediately. In case of swallowing rinse the mouth with water, drink a lot of water and consult a doctor.

Keep the person exposed at rest, do not force vomiting. Do not use in cases of known hypersensitivities to any of these ingredients.

Packaging

Portion of 30 g powder + 15 ml liquid

10260

For dental professional use. Made in Switzerland.

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© PD 200309