

# Thermosealer

Thermosealing machine

Equipment / Solutions for hygiene

Manual **thermosealing machine** for sterilization rolls.

## Characteristics

- High sealing speed
- Sealing area width 310 mm
- Intelligent validation thanks to the luminous warning indicator during use

## Advantages

- **Quality:** guarantees perfect sealing in a few seconds
- **Practicality:** the sealing lever is ergonomically designed and positioned frontally to simplify the packaging operations



Manual thermosealing machine



| Code      | Model                        |
|-----------|------------------------------|
| XR0150060 | Manual thermosealing machine |